Improved *in vivo* potency for divalent GalNAc-siRNA conjugates over tri-antennary GalNAc conjugates.

Adrien Weingärtner, Lucas Bethge, Marie Wikström Lindholm
Silence Therapeutics GmbH, Berlin, Germany
* m.lindholm@silence-therapeutics.com



1. Introduction

Targeting the hepatic asialo-glycoprotein receptor (ASGPR) with GalNAc-conjugated oligonucleotides has been shown to be sufficient for both, hepatocyte targeting and functional delivery as required for RNAi. Here we focus on the serial conjugation of serinol-GalNAc units to siRNA to investigate the impact of position and carbohydrate valency to achieve target knockdown *in vitro* and *in vivo*. In a systematic manner we probed dose response of siRNA conjugates at all four termini with GalNAc-units (n = 1 to 4) in primary murine hepatocytes. Moreover, we investigated distance dependencies for conjugates with two single GalNAc units. Finally, we assessed siRNA stability in endosomal and lysosomal milieu and saw clear improvements that correlate with improved *in vivo* duration and dose response.

9. Summary and conclusions

- \triangleright Simple design with two single serinol-GalNAc units positioned one at each terminus of the sense strand demonstrated 3-fold increased activity (ED₅₀ \approx 0.1 mg/kg) and duration of action compared to a tri-antennary GalNAc cluster, *in vivo*.
- Functional improvement by this design was confirmed in four out of four different sequences in vivo (data not shown).
- > This simple design results in increased RNA stability at low pH in endosomal & lysosomal environment.
- > Equal performance of two, three and four serial GalNAc units conjugated to siRNA, in vitro and in vivo.

2. Design and in vitro activity of serial GalNAc units at terminal siRNA positions A single GalNAc unit at either of the four termini of the siRNA showed strongly reduced or no dose response in silencing of target messenger RNA 24h post treatment in primary murine hepatocytes. Similarly, no silencing was observed for 5'-antisense conjugates with serial GalNAc units. However, at the 5'-sense, the 3'-sense and 3'-antisense siRNA conjugates with two and more serial GalNAc-units showed equal dose response of silencing as compared to the tri-antennary reference. ○ 2´OMe serinol-GalNAc tri-antennary



